

MRI Tables with fixed height

These only slightly magnetic stainless steel stretchers are used as patient transport tables in the field of magnetic resonance imaging (MRI, or MR). They are suitable for all MRI systems (according to their compatibility data sheets) up to **7 Tesla**.

The MRI patient transport tables can be used safely under the following conditions:

- Static magnetic field of 7 Tesla or less.
- With an additional spatial gradient of maximum 7 T/m or less,
- The table must not exceed the limit line of 7 T/m with the head end/foot end in the direction of movement towards the magnet.

Features

- The MRI transport tables with fixed height are available in standard versions with or without side safety rails.
- Four smooth-running locking castors (Ø 100 mm) for effortless patient transport.
- Continuously adjustable headrest (up to +40°) with clamping rod locking mechanism.
- Two-piece tabletop, upholstery thickness of 64 mm.
- · Maximum load 200 kg.
- Choice of fixed height 650, 700, 750 or 800 mm at no extra cost. The MRI tables will be supplied with a standard height of 650 mm, unless requested otherwise at time of ordering.
- Two standard widths available: 650 mm (Article No. U-3065/ LRVA and AWL-3065/VA) and 800 mm (Article No. EKG-3080 and AWL-3080/VA/LRVA).

Accessories

- Infusion pole, Article No. INF-1720/VA with four hooks and table fastening.
- Paper roll holder, Article No. P-3050/A (for table width 650 mm), or P-3070/A (800 mm) for mounting at the head end (not retrofittable).





Options

To ensure safe patient transport, the MRI tables are also available with side safety rails. In order to allow unobstructed access to the patient from either side as required, both side safety rails (800 x 250 mm, LxH) can be moved up or down by a pivot mechanism.



Article Number	U-3065/LRVA	EKG-3080/LRVA	AWL-3065/VA	AWL-3080/VA
Picture(s)	1	Similar 1	2, in 3	Similar 2
Width (mm)	650	800	650 (795 with side guards)	800 (975 with side guards)
Length (mm)	1900 (500/1400)			
Height (mm)	650/700/750/800			
Castors	Four locking castors, Ø 100 mm (Ø 125 mm also available, at extra charge)			
Side guards (mm)	- 800x250 (LxH), pivoting			