

## **Cardiology** Tilt Tables

The vein treatment tables have been redesigned specifically for the tilt table test, a medical examination method whereby the blood pressure of patients is evaluated according to changes in position. These cardiology tilt tables are available in two versions, **without echo-cutout** (Article No. KT-1062/E and KT-1060/EE) and **with echo-cutout** (Article No. EKAT-1062/E and EKAT-1060/EE).

The complete product descriptions with additional features, as well as the available options and accessories may be found in our current catalogue (issue 42) on pages 38, 40-41 and on our website www.agasan.com.



1 EKAT-1060/EE, LR/Z-50



2 KT-1060/EE, LR/Z-50, WM-30

AGA Sanitätsartikel GmbH Löhner Strasse 198b, 32584 Löhne, Germany Tel: +49 (0)5731 82283 · info@agasan.com

Article Number	KT- 1060/EE	KT- 1062/E	EKAT- 1060/EE	EKAT- 1062/E
	Without echo-cutout		With echo-cutout	
Picture(s)	In 2	3	In 1	Similar 1
Price (EUR)				
Frame	Lacquered			
Length (mm) Head/Main section	1900 500/1400			
Width (mm)	620 (680 with standard rails)			
Height adjustment (mm)	730-930 (Electric)	_	730-930 (Electric)	_
Fixed height (mm)	_	730	_	730
Echo-cutout (mm)	_		300x250 (LxW)	
Tabletop	2 sections			
Upholstery thickness (mm)	64			
Maximum load (kg)	150			
Tilt adjustment	-20°/+74°			
Headrest adjustment	+40° (Clamping rods)			
Paper roll holder	P-2065/A			
Accessories	Pages 40-41			



3 KT-1062/E

## www.agasan.com