

Therapy Tilt Tables / Standing Boards

These therapeutic tilt tables, or standing boards (Article No. KT-1080/EE and KT-10100/EE) can be used in various ways, for example for the treatment of patients with neurological diseases, head or back injuries, or paraplegic patients. They are ideally suited for patient mobilisation in neurology, surgery, physiotherapy and orthopedics.

The complete product descriptions with additional features, as well as the available options and accessories may be found in our current catalogue (issue 42) on pages 39, 40-41 and on our website www.agasan.com.



1 KT-1080/EE

- Upholstery cutout at the foot end and extra-large, removable footboard provide the patient with the necessary physical and mental security during treatment.
- Electric height and tilt adjustments by maintenance-free, extremely quiet electric motors controlled by air pulses.
- Height adjustment controlled by a foot switch. Drive specifications: 230V, 50/60Hz, 3,7A, 360W, IPX4, $v=14\text{mm/s}$, 7000N.
- Tilt adjustment controlled by a hand switch. Drive specifications: 230V, 50/60Hz, 1,8A, 360W, IPX4, $v=7\text{mm/s}$, 10000N.

Article Number	KT-1080/EE	KT-10100/EE
Picture(s)	1, 2	Similar 1, 2
Price (EUR)		
Frame	Lacquered	
Length (mm) Head/Main section	1900 500/1400	
Width (mm)	800 (860 with standard rails)	1000 (1060 with standard rails)
Height adjustment (mm)	520-960 (Electric)	
Tabletop	2 sections	
Upholstery thickness (mm)	64	67
Maximum load (kg)	130	
Tilt adjustment	0° / +90° (Electric)	
Headrest adjustment	+35° (Grid)	
Paper roll holder	P-2080/A	P-20100/A
Accessories	Pages 40-41	



2 KT-1080/EE