

Phlebology Tables with 2-pc. Tabletop AGA-PHLEBO-LIFT & AGA-PHLEBO-MAT

The three-section combined phlebology tables (presented in a different flyer) are available with a tabletop in two sections, for specific use in pure phlebology and dermatology.

The AGA-PHLEBO-LIFT (Article No. VEKO-1060/EE) with a tabletop in two sections features continuous electric height and tilt adjustment. The AGA-PHLEBO-MAT (Article No. VEKO-1062/3E) features continuous tilt adjustment by electric motors but has a fixed working height of 930 mm (in the horizontal position). The height adjustment cannot be retrofitted.

The maintenance-free and extremely quiet electric motors are controlled by air pulses via a foot switch. Drive specifications (each unit): 230V, 50/60Hz, 3,7A, 850W, IPX4, v=14mm/s, 7000N.



1 VEKO-1060/EE, LR/Z-50

The complete product descriptions with additional features, as well as the available options and accessories may be found in our current catalogue (issue 42) on pages 37, 40-41 and on our website www.agasan.com.

| | AGA-PHLEBO-LIFT | AGA-PHLEBO-MAT |
|---------------------------|----------------------------------|----------------|
| Article number | VEKO-1060/EE | VEKO-1062/E |
| Picture(s) | In 1 | 2 |
| Price (EUR) | | |
| Frame | Lacquered | |
| Length (mm) | 1900 | |
| Head/Main section | 500/1400 | |
| Width (mm) | 620 (680 with standard rails) | |
| Height adjustment (mm) | 730-930 (Electric) | — |
| Fixed height (mm) | — | 930 |
| Tabletop | 2 sections | |
| Upholstery thickness (mm) | 64 | |
| Maximum load (kg) | 150 | |
| Tilt adjustment | -20° / +74° | |
| Headrest adjustment | +40° (Clamping rods) | |
| Paper roll holder | P-2065/A | |
| Accessories | Pages 40-41 | |



2 VEKO-1062/E