

Phlebology Tables with 3-pc. Tabletop AGA-PHLEBO-LIFT & AGA-PHLEBO-MAT

The vein treatment tables AGA-PHLEBO-LIFT have been developed specifically for the requirements of phlebology. The goal is to support physicians and medical personal in providing efficient medical treatment thanks to the electric adjustments (in height and tilt), the specially-designed accessories and the ease of use.

The vein treatment tables are available with a two- or three-section upholstered tabletop. The three-section AGA-PHLEBO-LIFT (Article No. VEKO-1060/3EE) is also available without height adjustment. This model, AGA-PHLEBO-MAT (Article No. VEKO-1062/3E), has a fixed working height and only the tilting may be adjusted, using electric motors. The height adjustment cannot be retrofitted.

The vein treatment tables AGA-PHLEBO-LIFT and AGA-PHLEBO-MAT with a tabletop in three sections can be converted into gynaecology, proctology or rectoscopy examination chairs in just a few steps, by raising the backrest with a gas pressure spring and dropping the foot section using locking pins.



1 VEKO-1060/3EE

| | AGA-PHLEBO-LIFT | AGA-PHLEBO-MAT |
|---------------------------|-------------------------------|----------------|
| Article number | VEKO-1060/3EE | VEKO-1062/3E |
| Picture(s) | 1, 2 | 3 |
| Price (EUR) | | |
| Frame | Lacquered | |
| Length (mm) | 1900 | |
| Head/Seat/Foot section | 1020/365/465 | |
| Width (mm) | 620 (680 with standard rails) | |
| Height adjustment (mm) | 730-930 (Electric) | — |
| Fixed height (mm) | — | 930 |
| Tabletop | 3 sections | |
| Upholstery thickness (mm) | 64 | |
| Maximum load (kg) | 150 | |
| Tilt adjustment | -20° / +74° | |
| Backrest adjustment | +30° (Gas pressure spring) | |
| Leg section adjustment | Fold down (Locking pins) | |
| Paper roll holder | P-2065/A | |
| Accessories | Pages 40-41 | |



2 VEKO-1060/3EE



3 VEKO-1062/3E

The complete product descriptions with additional features, as well as the available options and accessories may be found in our current catalogue (issue 42) on pages 36, 40-41 and on our website www.agasan.com.